

Software History Signateq®

Version	Date	History
2.02	19.09.2025	- Fixes custom calibration compression/tension & tension entry calculation
2.01	14.02.2022	- Improved sensitivity adjustment capabilities for all three load directions - Support for new XENAX® ASCII-commands SQAC and SQCM (implemented in XENAX® firmware revisions 5.10B and higher)
2.00	10.12.2021	- Support for three calibration types "test report", "custom calibration" and "factory calibration"
1.04	30.08.2021	- Support for hardware revision 2.0 - General timing optimizations
1.02	26.08.2020	- First Release